

CORiTEC® 150i series

WHAT ARE THE DIFFERENCES?

	CORiTEC® 150i PRO	CORiTEC® 150i dry
Number of axes	5	5
Number of tool positions	10	10
Spindle power / rpm	860 W / 100.000 rpm	860 W / 100.000 rpm
Zero-Point-Clamping (manually)	✓	✓
Dry processing	✓	✓
Wet processing	✓	-
98 mm blank processing	✓	✓
CAD/CAM blocks / glass-ceramics	✓	✓ (only dry)
Premilled abutments	✓ (CoCr/Ti)	✓ (CoCr)
Metal processing (CoCr)	✓	✓
Control-PC	✓	✓
10" HD touchscreen	✓	✓
Integrated water tank & filter	✓	-
Chip tray	-	✓
Ionizer	✓	✓
LED statusbar	✓	✓
Auto-calibration set	optional	optional
Compressed air	3 bar external air required	3 bar external air required
Collet	3 mm	3 mm
Turning processing (universal fixture required)	✓	✓
CAM Software	iCAM V5 smart / iCAM HD (based on hyperDENT)	
98 mm blank holder	✓ included	✓ included
Open C-Clamp holder	optional	optional

	CORiTEC® 150i PRO	CORiTEC® 150i dry
Adapters		
Universal fixture for glass-ceramic and premilled abutments	available as accessory	available as accessory (only dry processing)
3-fold glass-ceramic holder (for 98 mm blank holder)	available as accessory	available as accessory (only dry processing)
6-fold glass-ceramic holder (for 98 mm blank holder)	available as accessory	available as accessory (only dry processing)
C-Clamp blank holder	available as accessory	available as accessory
Nt-Trading premilled blank holder 6-fold (for 98 mm blank holder)	available as accessory	available as accessory
Medentika premilled blank holder 6-fold (for 98 mm blank holder)	available as accessory	available as accessory
DESS premilled blank holder 6-fold (for 98 mm blank holder)	available as accessory	available as accessory
Dental Direkt premilled blank holder 6-fold (for 98 mm blank holder)	available as accessory	available as accessory
Zfx premilled blank holder 6-fold (for 98 mm blank holder)	available as accessory	available as accessory
Glass-ceramic adapter for universal fixture	available as accessory	available as accessory (only CoCr processing)
Nt-trading holder for universal fixture	available as accessory	available as accessory
Medentika holder for universal fixture	available as accessory	available as accessory
Dental Direkt holder for universal fixture	available as accessory	available as accessory
DESS holder for universal fixture	available as accessory	available as accessory
Zfx holder for universal fixture	available as accessory	available as accessory